Nephele systém

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
The screenshot shows a window titled "Nephele - Search in reports". It has three tabs: "Text retrieval" (selected), "Image retrieval", and "Search in data". Under "Text retrieval", there is a section "Search in text fields:" with six input fields: "Report no.:", "Author of report:", "Co-authors:", "Conservator's name:", "Artwork name:", and "Author's name:". Below this is a "Fulltext" section with a single input field. At the bottom of the search options is a "More options:" section with a label "Number of results in window" and a value of "10". At the very bottom are two buttons: "Search" and "Back".

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

The image displays two screenshots of the 'Nephele - Search in reports' application window. The left window shows the 'Image retrieval' tab, which includes two empty boxes for 'Image in white spectrum' and 'Image in UV spectrum'. Below these is a 'More options' section with a 'Number of results' input field set to 4. The right window shows the 'Search in data' tab, which features a 'Data' section with checkboxes for SEM/E, RTG-diff, FT-IR, and MALDI-T. Each checked option has associated input fields for chemical elements, material names, and peaks. Below the 'Data' section is a 'More options' section with 'Number of results in window' set to 10 and 'Allowance' set to 0. Both windows have 'Search' and 'Back' buttons at the bottom.

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
The image shows the 'Nephele - Insert new report' window. On the left is a tree view with the following structure: 'Heading' (expanded), 'General info', 'Analysed entities' (expanded, with 'Entity' selected), 'Methods of material research', 'Samples' (expanded, with 'List of samples' selected), 'Notes', 'Data' (expanded, with 'XRF', 'SEM/EDX', 'FT-IR', 'RTG-diffraction', and 'MALDI-TOF MS' listed), 'Files', 'Addendum' (expanded, with 'Images' and 'Texts' listed), 'Conclusion of material research', and 'User/Group rights'. Below the tree view is the 'alma' logo and a vertical stack of buttons: 'Add analysed entity', 'Add sample', 'Add image (addendum)', 'Add text (addendum)', and 'Remove selected node'.

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
Heading - Entity 1

Supplied information about artwork

Conservator's name:

Artwork name:

Date:

Author's name:

Technique:

Size:

Location/owner or artwork:

ID number:

Localization of samples extraction:

Artwork state, notes:

Photodocumentation of samples extraction

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

A screenshot of the 'Data - SEM/EDX' window. The 'New record' section contains input fields for 'Sample number', 'Layer number', 'Sample type' (with a dropdown menu showing 'Cross section', 'Fragment', 'Loose sample', 'Microtome slice', 'Thin section'), 'Analysis spec.' (with a dropdown menu showing 'Cross section', 'Fragment', 'Loose sample', 'Microtome slice', 'Thin section'), 'Artwork name', 'Coverage', 'Chemical element' (with 'Ag' selected), 'Wt %', and 'At %'. There are buttons for 'Add record', 'Clear table', 'Remove record', and 'OK'. Below this is a large table area labeled 'SEM/EDX' and a text field at the bottom labeled 'Measuring conditions and notes'.

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A screenshot of the 'Data - FT-IR' window. The 'New record' section contains input fields for 'Sample number', 'Layer number', 'Sample type' (with a dropdown menu showing 'Cross section', 'Fragment', 'Loose sample', 'Microtome slice', 'Thin section'), 'Analysis spec.' (with a dropdown menu showing 'Cross section', 'Fragment', 'Loose sample', 'Microtome slice', 'Thin section'), 'Artwork name', 'General association', 'Specified association', and 'wavenumber peak (cm⁻¹)'. There are buttons for 'Add record', 'Clear table', 'Remove record', and 'OK'. Below this is a large table area labeled 'FT-IR' and a text field at the bottom labeled 'Measuring conditions and notes'.

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The image shows two screenshots of dropdown menus. The left one shows the 'Material name' dropdown with 'kalcit' selected, and the 'Formula' dropdown with 'CaCO3' selected. The right one shows the 'Material name' dropdown with 'kaolinit' selected, and the 'Formula' dropdown with 'Al2Si2O8(OH)4' selected.

The screenshot shows the 'Data - RTG-diffraction' window. It has a 'New record' section with fields for 'Sample number', 'Layer number', 'Sample type' (set to 'Cross section'), and 'Analysis spec.' (set to 'Cross section'). There are also dropdowns for 'Artwork name', 'Material name' (set to 'kalcit'), 'PDF' (set to '5-0586'), and 'Formula' (set to 'CaCO3'). Buttons for 'Add record', 'Clear table', 'Remove record', and 'OK' are present. Below the form is a table header 'RTG-diffraction'.

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
The screenshot shows the 'Data - MALDI-TOF MS' window. It has a 'New record' section with fields for 'Sample number', 'Layer number', 'Sample type' (set to 'Cross section'), and 'Artwork name'. There is also a field for 'Analysis spec.' (set to 'vajíčko') and a dropdown for 'Material name'. Buttons for 'Add record', 'Clear table', 'Remove record', and 'OK' are present. Below the form is a table header 'MALDI-TOF MS'.

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





The screenshot shows a window titled "Nephele - User information". It contains the following fields and buttons:

- Login:
- Name:
- Surname:
- E-mail address:
- Section: Password change
 - Old password:
 - New password:
 - re-enter new password:
- Buttons: Send, Back

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The screenshot shows a window titled "Nephele - User Account Manage...". It contains the following buttons:

- New user account
- User list
- New user group
- User group list
- Back

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The screenshot shows the 'Nephele - User list' interface. On the left, a list of users includes 'admin' and 'guest'. The 'admin' user is selected. The main form contains fields for Login (filled with 'admin'), Name, Surname, E-mail, Password, and Re-enter password. Below these fields are buttons for 'Delete user', 'Reset', 'Send', and 'Back'. On the right, the 'User rights' section features a 'User group' dropdown menu set to 'admins' and a list of 13 checked permissions: Zobrazení majitele default..., Zobrazení inventárního čísla, Tisk záznamů, Zobrazení UV a SEM obr., Zobrazení interpretace vrstev, Zobrazení závěrů, Uzamčení záznamů, Fulltextové vyhledávání, Editace záznamů, Zobrazení obrázků pod lupou, Administrátor, Vložení nového záznamu, Přidávání nových látek pro data, and Mazání záznamů.

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